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
**OIL-FREE COMPRESSOR
COMPRESSEUR SANS HUILE
ÖLFREIE TROCKENLUFTKOMPRESSOREN
COMPRESOR DE AIRE SECO**

- MANUALE ISTRUZIONI
- OPERATOR'S HANDBOOK
- MANUEL D'UTILISATION
- GEBRAUCHSANWEISUNG
- MANUAL DE INSTRUCCIONES

CATTANI S.p.A.

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INDICE

Pagina

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INDEX

Page

— GENERAL RUNNING DATA	16
— RATED ELECTRICAL CHARACTERISTICS OF COMPRESSOR MOTORS	17
— LEGEND OF COMPONENTS	18
— INTRODUCTION	19
— SIGNS AND WARNINGS	19
— OIL-FREE COMPRESSOR AND COMPRESSED AIR DRYING SYSTEMS	19
— COMPRESSED AIR STERILIZATION	20
— LARGE PLANTS	20
— BREATHABLE AIR OR MEDICAL AIR	20
— FUNCTIONING OF COMPRESSOR	20
— CHOICE OF THE MODEL	21
— BLOK-JET	21
— ELECTRIC CONTROL PANELS AND TECHNOLOGICAL ROOM	21
— STARTING	22
— ORDINARY MAINTENANCE	23
— EXTRAORDINARY MAINTENANCE	24
— NOTICE	25
— TRANSPORT AND STORAGE	25
— ORDINARY MAINTENANCE RECORD	26
— EXTRAORDINARY MAINTENANCE RECORD	27
— TABLE OF COMPONENTS	70

General running data

Common features

Insulation class: I

Functioning: alternating service with physical-adsorption dryer

Operating temperature: from + 5 °C to + 30 °C

Motor rotating speed: at 50 Hz 1400 rpm - at 60 Hz 1600 rpm

Produced air at 5 real bar:

1-cylinder compressor 67.5 N l/min - 24 l tank

2-cylinder compressor 160 N l/min - 24/45 l tank

3-cylinder compressor 238 N l/min - 75 l tank

2-cylinder twin-head compressor 320 N l/min - 100 l tank

3-cylinder twin-head compressor 476 N l/min - 150 l tank

3-cylinder 3 head compressor 714 N l/min - 300 l tank

Blok-Jet:

6 heads 978 N l/min - 500 l tank

9 heads 1467 N l/min - 725 l tank

12 heads 1956 N l/min - 900 l tank

	Alternating current	IEC 417-5032
3N 	Three-phase alternating current with neutral	IEC 335-1
	Earthing (of functioning)	IEC 417-5019
PE	Protective conductor	CEI EN 60439-1
N	Neutral conductor	IEC 446
	Dangerous voltage	IEC 417-5036
	OFF (disconnection from mains)	IEC 417-5008
I	ON (connection to mains)	IEC 417-5007

The manufacturer is willing to provide suggestions, instructions, and to supply spare parts, literature, and any other useful information. The appliances are guaranteed for three years from date of sale, provided that guarantee card addressed to the manufacturer is returned to the same reporting date of sale, retailer's stamp and customer's name. Guarantee and manufacturer liability cease in case the appliances are found tampered by any kind of action performed by unfit and so unauthorised people.

Manufactured by ESAM S.P.A. - PARMA - ITALY

Rated electrical characteristics of compressor motors:

1-cylinder compressor

single-phase (1~) 50 Hz: 220 V - 0.55 kW - 3.7 A
 240 V - 0.55 kW - 3.7 A
 three-phase (3N ~) 50 Hz: 220/380 V - 0.55 kW - 3.1/1.8 A
 single-phase (1~) 60 Hz: 220 V - 0.65 kW - 4.6 A
 110 V - 0.65 kW - 9 A

2-cylinder compressor

single-phase (1~) 50 Hz: 220 V - 1.2 kW - 7.6 A
 240 V - 1.25 kW - 7 A
 three-phase (3N ~) 50 Hz: 220/380 V - 1.5 kW - 6/3.5 A
 single-phase (1~) 60 Hz: 220 V - 1.5 kW - 8.8 A
 110 V - 1.3 kW - 17.5 A

3-cylinder compressor

single-phase (1~) 50 Hz: 220 V - 1.5 kW - 10 A
 240 V - 1.5 kW - 9 A
 three-phase (3N ~) 50 Hz: 220/380 V - 1.5 - 6/3.5 A

2-cylinder twin-head compressor

single-phase (1~) 50 Hz: 220 V - two motors 1.2 kW - 7.6 A each
 240 V - two motors 1.25 kW - 7 A each
 three-phase (3N ~) 50 Hz: 220/380 V - two motors 1.5 kW - 6/3.5 A each
 single-phase (1~) 60 Hz: 220 V - two motors - 1.5 kW - 8.8 A each
 110 V - two motors - 1.3 kW - 17.5 A each

3-cylinder twin-head compressor

single-phase (1~) 50 Hz: 220 V - two motors 1.5 kW - 10 A each
 240 V - two motors 1.5 kW - 9 A each
 three-phase (3N ~) 50 Hz: 220/380 V - two motors 1.5 kW - 6/3.5 A each

3-cylinder 3-head compressor

three-phase (3N ~) 50 Hz: 220/380 V - three motors 1.5 kW - 6/3.5 A each

3-cylinder Blok-Jet

three-phase (3N ~) 50 Hz: 220/380 V - 6 motors 1.5 kW - 6/3.5 A each
 + dryer (1~) 50 Hz - 0.04 kW - 0.2 A
 three-phase (3N ~) 50 Hz: 220/380 V - 9 motors 1.5 kW - 6/3.5 A each
 + dryer (3 ~) 50 Hz - 0.1 kW - 0.4/0.2 A
 three-phase (3N ~) 50 Hz: 220/380 V - 12 motors 1.5 kW - 6/3.5 A each
 + dryer (3 ~) 50 Hz - 0.1 kW - 0.4/0.2 A

LEGEND OF COMPONENTS (table page 70)

A	1 st AIR FILTER (PAPER)
B	PISTON
C	CYLINDER
D	AIR-AIR EXCHANGER
E	DRYING COLUMN
F	AIR-DEW SEPARATING CYCLON
G	CYCLON WATER TANK
H	2 nd AIR FILTER (SINTERED BRONZE)
I	ADSORBENT SILICA GEL COMPOUND
L	3 rd FILTER (PAPER) - S STERILIZING FILTER
M	COMPRESSED AIR TANK
N	SELECTOR VALVE
O	SMALL AIR TANK FOR SILICA GEL REGENERATION
P	WATER COLLECTING BOTTLE
Q	ANTI - VIBRATION SUPPORTS
R	DISK FILTER
S	STERILIZING FILTER
T	MICROFIBER FILTER
1	ON SWITCH
2	OFF SWITCH
3	MOTOR
4	AIR-AIR EXCHANGER FAN
5	GAUGE INDICATING THE PRESSURE INSIDE THE TANK
6	MOTOR COOLING FAN
7	HEAD COOLING FAN
8	WATER CHECK TAP
9	DOOR
10	GRID
11	RED KNOB
12	DOOR LOCK RED SWITCH
13	SOLENOID VALVE

Introduction

The following presentation aims at illustrating the equipment and systems dealt with herein to users and engineers; **it also aims at explaining operation and maintenance**, as well as the dangers with the precautions required for accident prevention to users.

Signs and warnings

- **Danger of electric shocks: also 230 V ~ can be lethal.**


- **High temperature.**

- **General danger sign.**

- **Compulsory direction of flow or rotation.**

Signs cannot always fully express danger warnings, therefore it is necessary that the user reads the warnings and keeps them in due consideration. Failure to respect a danger sign or warning may harm operator or damage the equipment.

Do not remove protections, do not tamper with machines or their operation.

In particular do not carry out any welding or operations of any kind on the tank.

In spite of our best efforts, it is possible that the hazard warnings are not exhaustive; we apologise in advance to the user for this, and ask in the meantime that users anticipate any possible sources of danger that we might not have noticed, and to advise us of these.

Oil-free compressor and compressed air drying systems

When compressed air must be clean, hygienic, sterile or medical, the compressor must be a dry (oil-free) appliance and compressed air must be dried.

Water-oil emulsion, produced by a standard lubricated compressor, is known to be detrimental to the dentist's rotating instruments; moreover, since compressed air is used also to dry preparation before filling, even small traces of this emulsion can affect cementation. The matter is even more serious if compressed air is used in a sterile field: a sterilizing filter can effectively treat only dry air. Oil-lubricated compressors, fitted with filtering systems to retain humidity and emulsified oil, are more difficult to use, less safe and eventually more expensive than oil-free compressors.

Compressed air sterilization (draw. page 70)

Compressed air is filtered several times inside the circuit and conveyed through a silicagel compound which adsorbs humidity at a temperature within the range $+5\text{ }^{\circ}\text{C}/+35\text{ }^{\circ}\text{C}$. A pause is needed at every compressor charge in order to regenerate the silica gel (see below). Filtered and dried air can be considered suitable for instruments and hygienically safe for dentistry. Air quality can be improved with a sterilizing filter (S) to be installed at the compressor air outlet. The "SA" Balston filtering cartridge, with borosilicate glass microfiber filtering media (T), delivers a 99.9999+% efficiency on $0.01\text{ }\mu\text{m}$ particles. Balston filters can be sterilized in autoclave at $135\text{ }^{\circ}\text{C}$ (max. 20 cycles); when used with clean air, life of the filters is about 1 year, unless operating conditions require more frequent substitutions. At this stage the air can be considered as pure, if conveyed into a reasonably short piping which is not subject to thermal stress.

Large plants (draw. page 74)

Longer piping, such as those you find in hospitals, clinics, or schools, require a different air drying system. The recommended air treatment is that of selective membrane, or silica gel double column with the automatic regeneration and alternation of columns; these systems guarantee a proper drying of the air, but do not require pauses for regeneration. For large plants we can supply various versions of Blok-Jet units (see pag. 74), with membrane dryers (produced air: 978-1467-1956 N l/min).

Breathable air or medical air

The oil-free compressor is ideal for producing breathable air for hospitals. In fact, in absence of oil, no gas that may harm patients is released, therefore evidently the hazardous gas filtering system releasing the air produced by the compressor will be less affected. **To obtain medical air after the compressor, air goes through a series of filters removing all chemical and physical impurities, in compliance with the provisions of the monograph Issue X of the Italian Official Pharmacopoeia.** More specific and detailed information on medical air is available upon request.

Functioning of compressor (draw. page 70)

The air enters the cylinder through filter (A) and is filtered for the first time; this filter traps atmospheric dust, ensuring air cleanliness and thus protecting pistons and cylinders. Piston (B) slides inside the cylinder (C) without any lubrication; material coating the piston is compatible with the liner, and in spite of friction both liner and piston are not subject to any significant wear.

The air is compressed and warmed in the cylinder, therefore it is necessary to convey it through an air-air exchanger (D) where, at a room temperature of about $20\text{ }^{\circ}\text{C}$, ΔT is stable at about $5\text{ }^{\circ}\text{C}$. Air temperature needs to be reduced to dew point in order to collect condensation.

Cooled air is then conveyed through cyclon (F) in the drying column (E) thus forming dew, which is collected in the tank (G). At the top of the column, the air goes through a disk filter (R), and then through a high-adsorption silica gel (I) compound and is completely dried. Two filters are located at the drying column air outlet: the first one is in sintered bronze (H), the second one in paper (L). Upon request, we can supply a Balston sterilizing filter (efficiency: 99.9999 +% for $0.01\text{ }\mu\text{m}$ particles). Then the compressed air, dried and hygienically filtered, enters the tank (M) to be used. The tank is coated with certified

alimentary resin guaranteeing good conservation of air. Whilst the tank is being filled, the air also enters a small reservoir (O) through the selector valve (N). At the end of each cycle, the dry air in the small reservoir flows back through the drying column in the opposite direction and regenerates the silica gel compound. This reversed flow of air carries all the water (previously removed from the compressed air) and drains it into the bottle (P). The whole process is automatic.

Choice of the model (draw. page 70-71-72-73)

The type of oil-free compressor will be selected according to the dentist's rotating instruments and to the compressed air demand of the surgery, clinic, hospital or factory. One cylinder of the compressor produces from 60 to 80 l/min, approximately, while the dentist instruments subject to continuous use generally do not consume more than 60 l/min of air at the pressure of 5 bar, so one cylinder can adequately supply one chair. The capacity of the compressor must be abundant, in any case, with respect to consumption, to facilitate the pauses, necessary to regenerate the drying column at the end of the charging operation. The double column set with automatic drainage does not need any pause. The loss load that occurs along the distribution line is generally due to the lack of a ring type layout and to pipes with an excessively small cross section. All too often these important aspects are neglected. As for the noise level, the compressor can be supplied with plastic cover or deadening box (see pag. 71-72-73), allowing for a noise reduction of 10 and 20 dB (A) respectively.

The pressure of the one-cylinder-compressor tank is set by the manufacturer at min. 5.5 bar and max. 7.5 bar. Regulation of two-or three-cylinder-compressors with ranges from 6 to 8 bar. Different regulations of the working pressure will affect the life of the compressor: the higher the max pressure the greater the effort and wear of the machine. It is prudent, in any case, to keep the maximum pressure at least 1.5 bar below the level of pressure of the safety valve, set at 10.5 bar.

Blok-Jet (draw. page 74)

Blok-Jet consists of an enclosure housing aspiration and compression systems, or only a series of compressor heads. The enclosure allows to place one appliance on top of another (compressors are usually placed on top of aspirators), thus saving 50% space. We can supply customized Blok-Jet units, once surgery or clinic, hospital or factory requirements and available space are known. According to location and customer demand, Blok-Jet can be produced in different versions: open, side-panelled, fully panelled and ventilated, and also with soundproof conditioned box. A group of heads, controlled by an electronic board which starts the heads in sequence, is more versatile and less expensive than a single appliance, particularly for large plants. Groups of 6-9-12 heads and more, (draw. page 74) allow replacing one head without switching off the plant. Sequential starting of heads allows to avoid dangerous peak amperage.

Electric control panels and technological room

It is advisable to install suction units, compressors, electric panels controlling machine operations, amalgam separators with possible filtering lines and drain pumps, in a duly equipped technological room. Check that the air exhausted from the suction unit is conveyed outside and that it is not sucked up by the compressor.

Keep the machine room clear of anything that is not related to the machines contained therein, taking special attention of flammable materials; make sure that there is no

risk of the formation of corrosive, flammable or explosive mixes. Install a smoke detector for fire prevention, connected to a non-stop surveillance service.

Make sure that in this room the temperature is a minimum of +5 °C and a maximum of +35 °C. Electric supply must be sufficient for the loads specified on machine plates and three phases + neutral + ground must be in place. Grounding is an indispensable condition also for single-phase supply. Supply lines, electric units and machines must be protected from direct and indirect electric contacts, overloads and overcurrents, in compliance with C.E.I. 64-8 (IEC correspondents available on request) regulations for first class equipment or, in compliance with C.E.I. 64-4 (electrical plants in medical rooms). It is advisable to protect the supply network from extra-currents. Protect electric panels and machines from accidental sprinkles, provide for the collection of liquids on the ground with drainage connected to the sewage system. Outside the room, install a fire extinguisher of appropriate size and periodically check its efficiency. Allow access to personnel only, make sure that all the personnel wear appropriate protections from sprinkles and accidental contacts with infected parts; ensure that no repair or inspection is carried out on machines in operation or ones connected to the mains. Machine boxes are not only a means to reduce sound vibrations, but it is also a protection from accidental contacts and (infected) sprinkles from machines under pressure.

Employ boxes produced by the compressor manufacturer only. The manufacturer cannot be liable, even within the warranty period, for units fitted with inappropriate boxes.

The machine installation must be carried out by an expert operator, duly trained and supplied with original spare parts. The operator assembling the equipment must refer to the machine manual to carry out the machine final test and to instruct users (surgery staff) on use and routine maintenance.

Technical maintenance must be carried out by engineers duly trained and supplied with original spare parts.

Arrange for a periodical inspection of equipment and systems. This inspection is not only a means to prevent the surgery from forced stops and to check its operating conditions, but it is also a way of preventing accidents and fires (According to Italian law 626).

Starting (draw. page 70)

- **Danger of electric shocks: also 230 V ~ can be lethal.**

- **High temperature.**

- **General danger sign.**

- **Compulsory direction of flow or rotation.**

- **Unpack the compressor following the instructions shown on the package.**
- **Dispose of the package in compliance with regulations.**
- **Verify that the appliance has not been damaged during transport.**
- **Do not connect damaged appliances to the mains.**
- **Do not use extension leads, multiple plugs or sockets.**
- **Ascertain that the feeding line is adequate to feed the compressor.**
- **Ascertain that the compressor (or Blok-Jet) is installed in a clean spot, far from heat sources and from depots of contaminating substances, and that it aspirates clean air, free from dust, gases, and humidity.**

Before starting the compressor, ascertain that piping is clean; heavy debris could damage the appliances.

When the compressor is connected to the mains and to the distribution line, push on the black button (1) to start it. In case of a three-phase compressor, check the direction of rotation marked on the motor (3) and on the grid (10). Check that the ventilator (4) of the air/air exchanger works properly.

As soon as the manometer (5) reads 8 bar*, the compressor stops and will start again at 6 bar. At every stop a slight whistling warns that the silica-gel is being regenerated. Fans (6 and 7), installed on the same shaft, turn simultaneously with the motor. As a rule, the hottest points of a properly running and pausing compressor get temperatures ranging from 60 to 85 °C.

In one, two or three cylinder compressors, in case of any anomaly, press the red knob on yellow background (11) placed on the door (9) (this door should be kept closed when the compressor is running) or the red button (2), (if the door had been left unintentionally open). For tandem compressors or Blok-Jet, turn the red on yellow door-locking switch (12) placed on the control panel. For special plants, clinics, hospitals and factories, customized manuals are issued.

Ordinary maintenance (drw. page 70)

To be entrusted to the surgery staff, purposely trained

Before approaching the compressor, disconnect it from the mains and ascertain that nobody can connect it again without the maintenance operator knowing it.

If the compressor is installed in a dusty spot or where it can aspirate solid particles, as dust, sand, dry leaves or similar, the head filters should be cleaned and replaced very often.

If it is installed in a clean room, it will be enough to clean the filters (A) every three months and replace them every six months. If a Balston sterilising filter is installed (S), we recommend to sterilise it in autoclave at 135 °C ∞ at least every six months, for max 20 times. We recommend to wear disposable gloves and to put the filter immediately into the autoclave.

A Balston filter lasts one-year max. Replacement cartridges are available; follow the directions given with the filter and mark the operation on the report sheet.

Check the ventilator running every week, as the compressor cannot work for long without adequate ventilation.

Tap (8), placed under the air tank, must be opened at regular intervals: if moisture is noticed, call for an engineer. Every week, loose manually the condensation draining screw placed under the tank: air coming out must be dry; otherwise call for an engineer.

Note that the air-drying system cannot work properly with environment temperature above + 35 °C.

The compressor should not start if no compressed air is used: differently, check for any leak in the system or in the compressed-air-using appliances.

Charging time is about 45/55 sec. and the pausing for the regeneration of the silica gel and for cooling time is about one third of the charging one.

A careful watch lengthens the life of the appliance: when it gets noisy or vibrates, something has loosened or worn out – call for the engineer.

Experience and the volume of the surgery work will give indications to every operator about increasing or decreasing the frequency of the mentioned operation in comparison with our advice.

Always fill in the “Ordinary maintenance” report-sheet.

* Excepted one-cylinder compressor running from 5.5 to 7.5 bar.

Extraordinary maintenance

To be entrusted to a trained and authorized engineer, provided with original spare parts.

Checks must be regular; their frequency is related to the volume of the surgery work. A surgery working 8 hours a day, for five days a week, needs one check every three months, if appliances are watched also by the surgery staff charged with ordinary maintenance. The engineer charged with the extraordinary maintenance must use original spare parts only, must not modify the appliances or their functioning, and must not modify any safety device. In particular he must not carry out any welding on the compressor tank.

Before any operation, consult the instruction manual, split-up drawings, and electrical diagrams.

Before approaching the compressor, disconnect it from the mains; if the main switch is far away and cannot be surveyed closely, lock it. On the first check, after some functioning of the compressor, the engineer must check that no terminal of the control panel has loosened and that no overheating signs are visible on terminals and cable connectors. He must also check that the cylinder screws are not loosened - use a dynamometric key adjusted at 6 N•m. Ascertain that the compressor gets the set max. pressure at every charge. Charging time is about 45/55 sec. and the silica gel regeneration and cooling time is about one third of the charging one.

If the compressor has difficulty in getting the max set pressure, check the tension on the line and capacity of condenser. Check the electrical absorption.

At every stop of the compressor a dry-air-reversed-flow dries the silica gel; if it is not so, check or replace the electrovalve (13). If the compressor has difficulty in getting the working rate, check the tension on the line and capacity of condenser. Every six months check the absorption with ammeter pliers. Check for any leak on the appliance, on the line, and on the dental units.

Check the aspiration value at every cylinder (from 6 to 8 bar): aspirated air must not be less than 6000 N l/h. Replace burnt bulbs in the control panel, replace flaming relays and remote switches and those with worn-out contacts. Do not modify the functioning or the electrical and mechanical protections. Any change of the running noise can be a sign of malfunctioning and of a breakdown risk, therefore it is a good rule to replace noisy heads.

Check the engine room temperature: the air-drying system cannot work properly at temperatures above + 35 °C.

Always fill in the “Extraordinary maintenance” report-sheet.

Notice

- The manufacturer is willing to supply spare parts, technical information and any other useful information.
- Distributors, agents retailers, and authorized service engineers are supplied with split-up drawings, electrical diagrams, handbooks and updating, as regards servicing and maintenance.
- The appliance is guaranteed for 3 years from date of sale, provided that guarantee card addressed to the manufacturer is returned to the same reporting date of sale, retailer stamp and customer's name.
- Guarantee and manufacturer liability cease in case appliances and/or plants are found tampered by any kind of action performed by unable and thus unauthorised people.
- For any use not contemplated or specified in this handbook please refer to the manufacturer.
- On the web-site: www.cattani.it you can find our **updated manuals**. We suggest to consult them especially concerning the **security**.

Transport and storage

- *Packed equipment can be transported and stored at a temperature range of - 10 °C + 60 °C.*
- *Packages must be kept away from water and splashing and cannot tolerate humidity >70%.*
- *Packages with the same weight can be stored in piles of three only.*

ORDINARY MAINTENANCE RECORD

Compressed air distribution system

Record to be filled in at every ordinary maintenance operation.

The ordinary maintenance must be entrusted to the surgery staff purposely trained.

Photocopy this record and keep it after filling in.

Date of maintenance Operator's signature		
Before getting near the compressor, cut the power off	YES	NO
Every three months: clean the filters on the head (in a clean air environment)*	YES	NO
Weekly: check the functioning of ventilators	YES	NO
Weekly: check the presence of a coming out air blow at every stop	YES	NO
Monthly: Check the temperature in the technical room	YES	NO
Monthly: loose the screw at the bottom of the tank, ascertain that the coming out air is dry, and turn again the screw	YES	NO
Every six months: replace the filters on the head (in a clean air environment)*	YES	NO
Keep the technical room clear of anything that is not related to the machines contained therein, taking special attention of flammable materials; make sure that there is no risk of formation of corrosive, flammable or explosive mixes	YES	NO
At every maintenance operation: check there is no hindrance (dust or other) to the proper ventilation of the machines	YES	NO

** In dusty environments clean and replace the filters more often, according to environmental conditions*

EXTRAORDINARY MAINTENANCE RECORD

Compressed air distribution system

Record to be filled in at every extraordinary maintenance operation.
The extraordinary maintenance must be entrusted to a trained engineer, provided with original spare parts .

Photocopy this record and keep it after filling in.

Date of maintenance Operator's signature		
Before getting near the compressor, cut the power off	YES	NO
At the first check up: check the screws of the terminals in the control panel	YES	NO
At the first check up: verify the head screw with a dynamometric key 6N.m	YES	NO
Every three months: clean the filters on the head (in a clean air environment)*	YES	NO
Every three months: verify the charging time (from 6 to 8 bar = 45/55 sec.)	YES	NO
Every three months: check the system and seek for any leak	YES	NO
Every three months: sterilize the absolute filter by autoclave	YES	NO
Every three months: verify the sound and visual signals	YES	NO
Every six months: replace the filters on the head (in a clean air environment)*	YES	NO
Every year: replace the absolute filter on the column	YES	NO
At every check-up: verify the functioning of ventilators	YES	NO
At every check-up: verify the functioning of drying system	YES	NO
At every check-up: verify the temperature of the technical room	YES	NO
At every check-up: check the noise	YES	NO
At every check-up: check the capacity of the capacitor	YES	NO

** In dusty environments clean and replace the filters more often, according to environmental conditions*

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











GAMMA COMPRESSORI

COMPRESSOR RANGE



MONOCILINDRICO
ONE-CYLINDER COMPRESSOR



BICILINDRICO 25L
TWO-CYLINDER COMPRESSOR 25-litres tank

BICILINDRICO 50L
TWO-CYLINDER COMPRESSOR 50-litres tank



TRICILINDRICO
THREE-CYLINDER COMPRESSOR



GAMMA COMPRESSORI

COMPRESSOR RANGE


BICILINDRICO TANDEM
TWIN-CYLINDER TWIN-HEAD COMPRESSOR





TRICILINDRICO TANDEM
THREE-CYLINDER TWIN-HEAD COMPRESSOR

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ITALIAN PATENTS OR PATENT APPLICATIONS:

CATTANI: 1234828 - 1259318 - 1.187.187 - 1253460 - 233634 - 2337706 - 1294904 - 1305443 - 1310808 - 1340812
ESAM: 1225173 - 1253783 - 0791751

FOREIGN PATENTS OR PATENT APPLICATIONS:

CATTANI: US 4,787,846 - US 5,039,405 - US 5,002,486 - US 4,684,345 - US 5,330,641 - CH 0211808 - DE 0211808
- FR 0211808 - GB 0211808 - SE 0211808 - DE 0335061 - ES 0335061 - FR 0335061 - GB 0335061 - AT 0557251
- DE 0557251 - ES 0557251 - FR 0557251 - GB 0557251 - DE 0638295 - DK 0638295 - ES 0638295 - FR 0638295
- GB 0638295 - NL 0638295 - SE 0638295 - US 6,083,306 - US 6,090,286 - US 6,022,216 - US 6,423,124
ESAM: US 4,948,334 - DE 0351372 - ES 0351372 - FR 0351372 - GB 0351372 - US 5,779,443 - CH 0791751 - DE 0791751
- ES 0791751 - FR 0791751 - GB 0791751 - PT 0791751 - ES 107358 - FR 222.394/395

PENDING PATENT

CATTANI: EP 99830010.7 - EP 99830011.5 - EP 99830250.9 - EP 00830491.7 - EP 02425482.3
ESAM: EP 02425317.1 - MO 2001A000176 - US 2003-0044295

CATTANI s.p.a.

VIA NATTA, 6/A - 43100 PARMA - ITALY

TEL: +39 0521 607604 - SALE DEPT. FAX: +39 0521 607628

PURCHASING DEPT. FAX: +39 0521 607855 - ACCOUNTING DEPT. FAX: +39 0521 399966

<http://www.cattani.it> Email: cattani@tin.it

Company with Quality System Certified by DNV UNI EN ISO 9001/2000 - UNI CEI EN ISO 13485

esam s.p.a.

VIA NATTA, 4/A - 43100 PARMA - ITALY

TEL: +39 0521 607613 - SALE DEPT. FAX: +39 0521 607628

PURCHASING DEPT. FAX: +39 0521 607855 - ACCOUNTING DEPT. FAX: +39 0521 399966

<http://www.esam.it> Email: esamsa@tin.it

Company with Quality System Certified by DNV UNI EN ISO 9001/2000